

Mobile Designs, Inc.

A PROFESSIONAL AND MODERN LOOK IN DONOR BED

MD2500-3 MOBILE BLOOD DONOR BED (3 beds with transport carrier)



Specifications:

| Description | Height | Width | Depth | Weight |
|-----------------------------|---------------|--------------|----------------|----------------|
| 1 Lounge, Folded | 109 cm | 60 cm | 26 cm | 12.7 kg |
| 1 Lounge, Opened | 136 cm | 60 cm | 130 cm | 12.7 kg |
| 3 Lounges On Carrier | 115 cm | 60 cm | 91.5 cm | 47 kg |

Standard features:

- ❖ Comfortable and foldable bed for donor laying on during mobile blood donation
- ❖ 4 position mobile bed with positive lock in the lay back positions
- ❖ Aluminum frame construction
- ❖ Arm rest swings out of the way for easy access during blood donation
- ❖ Arm rest also adjusts up and down and can be used on either side of the bed
- ❖ Seat & back covers material made of are fire retardant & disinfectant resistant textaline
- ❖ All covers are removable for cleaning, repair or replacement
- ❖ Head rest is adjustable for the individual
- ❖ Mobile bed load capacity – 181.5 kg
- ❖ Each unit of bed includes: 1 arm rest, 1 head rest and seat cushion
- ❖ Transport carrier able to accommodate three units of bed during transportation or storage

Optional accessories:

- ❖ Second arm rest (Part No. MD-2501)
- ❖ Scale mounts (Part No. MD-2502)
- ❖ Back pouch (Part No. MD-2505)
- ❖ Arm rest pouch (Part No. MD-2503)
- ❖ Transport carrier (Part No. MD-2504)
- ❖ Accessories pouch, for 6 arm rests (Part No. MD-2506)

Manufactured by

MOBILE DESIGNS, INC.

4650 CATERPILLAR ROAD

REDDING, CA 96003

PHONE: 530-244-1050 FAX: 530-244-5016