

---

# *Mobile Designs, Inc.*

A PROFESSIONAL AND MODERN LOOK IN DONOR LOUNGES

---

## MD 3500P DONOR LOUNGE

(with Pump)



### **SUPPORT AND COMFORT**

This reclining lounge can be tilted forward for easy entry. It can also be tilted further back than most other lounges making it a lot easier to treat donor reactions. The lounge includes an adjustable headrest that supports your upper shoulders, neck and head. The armrest is heavily foamed so your elbow sinks into the cushion for more comfort. It adjusts in or out from the chair, up & down and swings out of the way for easy access. This lounge is built with quality materials to provide long life and dependable service. Each lounge includes 2 armrests and 1 headrest.

### **SPECIFICATIONS:**

Seat Height	: 61 – 76 cm (Adjustable)
Arm Rest Height	: 76 – 102 cm (Adjustable)
Scale Height	: 46 – 61 cm (Adjustable)
Height	: 104 – 145 cm (Adjustable)
Width	: 89 cm (Inside armrests)
Depth	: 188 cm
Capacity:	: 180 kgs

### **OPTIONAL ACCESSORIES:**

MD 3504	DVD Holder
MD 3502	Scale Mount

---

Marketed by:

***Biomarketing Services (M) Sdn Bhd.***

No 21, Jalan 4/62A, Bandar Manjalara, Kepong, Kuala Lumpur 52200

Phone: +603-6273-3068 Fax: +603-62720093

Email : [info@biomarketing.com.my](mailto:info@biomarketing.com.my)