

TECHNICAL DATA

MODEL	COMFY COUCH C102
Input power	220 ± 10% ~ VAC, 50 Hz, Single Phase
Power	120 W
Classification	Protection against electrical shock, Class I, Type B.
Output Socket Voltage	220±10% V, 50Hz
Output Socket Power	60W Max
Fuse	250V, 1A, 5X20 Type M. Make: Little fuse.
Movement / Position	Any position from seating to Vasovagal attack position. Height adjustment from 620mm to 820mm.
Movement Control	By logically placed buttons on handheld wired remote.
Movement Actuation	By DC motor and Torque Multiplier
Maximum Load	150 kg.
Seat Width	520mm
Profile length	1920mm
Seat height	620mm
Width of arm rest	150mm
Length of arm rest	480mm
Wheel diameter	100mm
Weight of Equipment	90kg
Upholstery	4" soft foam with Rexene covering.
Accessories	
BCM tray size	400mm X 280mm
Sphygmomanometer tray size	270mm X 130mm
Accessory tray size	370mm X 160mm
IV Stand	Height adjustable IV Stand.